



METHOD-1 SHIELD™

ANTIMICROBIAL PROTECTION

PROTECTING THOSE WE COUNT ON

product of **Expo EXPO Displays®**

Using our Method-1 Shield, an antimicrobial laminate, we can now install protective graphics to your exhibit surfaces to help protect the people who touch them.



KIOSKS



RAILINGS



SEATING



COUNTERS & PODIUMS



TABLES



DOOR HANDLES


WHY USE METHOD-1 SHIELD?

Method-1 Shield graphics offer branding opportunities with effective and long-lasting protection against many forms of bacteria.

It's peace of mind for you, your staff and your trade show attendees.



Protect your trade show attendees and employees by protecting the areas they touch.



Get added exposure for your company in places you had not previously thought about.



Show your trade show attendees and employees you are doing everything possible to keep them healthy

PROMOTING YOUR PROTECTED SURFACES



FREESTANDING SIGNAGE

Make sure your trade show attendees realize that your surfaces are protecting them! We offer a customized 3.75" x 4.75" sign in a sturdy aluminum frame. This is the perfect way to offer reassurance without saying a word.



SURFACE DECAL

Look for this decal to identify where Method-1 Shield is applied. This easily recognizable label will provide visitors peace of mind knowing that your surfaces are safe to touch!



WINDOW/DOOR DECAL

Easily recognizable on the doors and windows of your exhibit conference room, this decal informs everyone of your advanced level of protection. This decal is available in 4"x4", 5"x5" and 6"x6" sizes.

METHOD-1 SHIELD GRAPHICS

Method-1 Shield is not just a protectant, it is also a graphic! So your protected surface can be any image, color or design imaginable. We can also add your logo to the signage and decals promoting the protection for additional branding for your company.



PROFESSIONAL INSTALLATION

ExpoDisplays recommends and offers professional installation of Method-1 Shield by our trained installation team. We can install Method-1 Shield prior to shipping your exhibit. You never have to give it a second thought.



**METHOD-1
SHIELD™**

ANTIMICROBIAL PROTECTION

PROTECTING THOSE WE COUNT ON

product of **Expo** **EX** **PO** **Displays**®

WHAT IS METHOD-1 SHIELD?

Method-1 Shield is a technically sophisticated and high-quality film, which helps provide effective and long-lasting protection for your public spaces, buildings, and public transportation. It can be applied under any circumstances and in all locations with high sanitary requirements, such as such as clean rooms, cold rooms, nurseries, hospitals, food industries, and trade show venues. Method-1 Shield is constructed out of PVC, which is then enriched with Silver Ions.

PROOF FROM TESTS: THE ACTIVITY OF METHOD-1 SHIELD

STRAIN	REDUCTION IN BACTERIA AS A %	LOGARITHMIC DECREASE
Salmonella enterica subsp enterica	99,99	> 4,6
Listeria monocytogenes	99,99	> 4,2
Staphylococcus aureus	99,99	> 4,1
Escherichia coli	99,99	> 4,5
MRSA	99,99	> 3,5
P.aeruginosa	99,99	> 5,7

Silver Ions are highly reactive and destroy multiple sites in bacteria cells. During cleaning operations, as well as in the presence of humidity, silver ions are released from the top layer of the film. These ions come into contact with the bacteria, blocking their metabolism and/or interrupting their proliferation mechanism, leading to their destruction.



Inhibits the development of 99.99% of germs tested.



Waterproof, smooth, and easy to wash material.



Effective even at low concentrations, perfectly ecological



Resistant to most chemical agents: alcohol, bleach and germicidal wipes.



Prolonged protection, 24 hours a day, 7 days a week, for several years.



Non-irritable and safe against human skin.

For for technical information and specifications visit:

[Expodisplays.com/method-1-shield/](https://expodisplays.com/method-1-shield/)

WHERE CAN METHOD-1 SHIELD BE APPLIED?

ANYWHERE WITH A SMOOTH SURFACE!

Method-1 Shield can be applied anywhere and everywhere you'd like to protect, even the most high traffic areas of your exhibit! Method-1 Shield is waterproof and the adhesion is immediate. After about 24 hours of contact the adhesion becomes permanent. Method-1 Shield is already used in a variety of locations including Trains, Planes, and Subway systems; thus exposing the film to a significant amount of constraints it is built to handle.

THE POTENTIAL USES ARE ENDLESS

There are many high touch areas within your exhibit that rarely get cleaned during the course of a regular trade show day. These areas are often cleaned only once per day. Protocols will most likely now call for multiple cleanings throughout the day in these publicly accessed areas.

Covering these type areas and items with Method-1 Shield will increase consumer confidence in the cleanliness of the areas and items that are touched most often. Seeing the Method-1 Shield label offers reassurance that everything that can be done is being done to keep everyone safe and healthy.



FOOD AND BEVERAGE STATIONS



INTERACTIVE TOUCH SCREENS WITHIN YOUR EXHIBIT



PRODUCT DISPLAYS WITHIN YOUR EXHIBIT